Instructions for Use

STRATAFIX™ Spiral MONOCRYL™ Plus

Knottless Tissue Control Device

Antibacterial

Undyed Monofilament Synthetic

Absorbable Unidirectional Tissue

Description

STRATAFIX™ Spiral MONOCRYL™ Plus Knottless Tissue Control Device is a sterile, antibacterial monofilament, synthetic absorbable device composed of 100 weight percent copolymer of glycolide and ε-caprolactone. The empirical formula of the copolymer is (C₄H₆O₂)_x(C₈H₁₂O₃)_y. Poliglecaprone 25 copolymer has been found to be nonpyrogenic and elicits only a slight tissue reaction during absorption.

STRATAFIX™ Spiral MONOCRYL™ Plus Device contains IRGACARE™ MP (triclosan), a broad spectrum antibacterial agent, at no more than 2360 µg/m. STRATAFIX™ Spiral MONOCRYL™ Plus Device, Variable Loop Design consists of barbed suture material, armed with a surgical needle on one end and a fixation loop at the opposite end. STRATAFIX™ Spiral MONOCRYL™ Plus Device is designed to anchor with a closed loop at one end and a unidirectional barbed section on the other end. STRATAFIX™ Spiral MONOCRYL™ Plus Device bars are oriented in one direction to allow tissue approximation without the need to tie surgical knots.

While the formation of bars in STRATAFIX™ Spiral MONOCRYL™ Plus Device reduces the tensile strength relative to non-barbed suture material of the same size, tying of knots in non-barbed suture materials also reduces their effective strengths. For this reason, the strength of STRATAFIX™ Spiral MONOCRYL™ Plus Device can be compared to United States Pharmacopeia (USP) knot strength of non-barbed sutures. USP designations for diameter are used to describe STRATAFIX™ Spiral MONOCRYL™ Plus Device suture material after barbing, except for minor variation in suture diameter with a maximum average of 0.1 mm.

STRATAFIX™ Spiral MONOCRYL™ Plus Device is available in a range of gauge sizes and lengths, attached to needles of various types and sizes, and in presentations as described in the How Supplied section.

The European Pharmacopoeia (Ph. Eur.) recognizes units of measure Metric and Ph. Eur. sizes as equivalent, which is reflected on the labelling. The USP and Ph. Eur. sizes of STRATAFIX™ Spiral MONOCRYL™ Plus Device are further defined in Table 1.

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STRATAFIX™ Spiral MONOCRYL™ Plus Device is intended for use only by healthcare professionals who are trained in surgical suturing techniques including the use of STRATAFIX™ Spiral MONOCRYL™ Plus Device.

Clinical benefits expected from soft tissue approximation are to promote wound healing by primary intention.

A summary of safety and clinical performance can be found at the following link (upon activation) : <https://ec.europa.eu/tools/edumed>.

Tensile Strength

STRATAFIX™ Spiral MONOCRYL™ Plus Device straight tensile strength is stronger than the knot tensile strength for a USP and EU Pharmacopeia poliglecaprone 25 device of the equivalent size as shown in Table 2.

Table 2. Tensile Strength Comparison

USP Device Size Designation	Ph. Eur. Device Size (Metric/Ph. Eur.) Designation	Device Minimum Knot Tensile Strength	Metric/Ph. Eur. (N)
2-0	3	2,68	26,8
3-0	2	1,77	17,5
4-0	1,5	0,95	9,5

Indications/ Intended Use

STRATAFIX™ Spiral MONOCRYL™ Plus Device is indicated for use in soft tissue approximation where the use of absorbable suture is appropriate.

Contraindications

STRATAFIX™ Spiral MONOCRYL™ Plus Device is not to be used where extended approximation of tissue under stress is required such as in fascia closure. This device is not to be used in conjunction with or for fixation of prosthetic devices (e.g., heart valves or synthetic grafts) that are non-absorbable in nature. STRATAFIX™ Spiral MONOCRYL™ Plus Device should not be used in

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patients with known allergic reactions (local hypersensitivity) to IRGACARE™ MP (triclosan).

Warnings

The safety and effectiveness of STRATAFIX™ Spiral MONOCRYL™ Plus Device have not been established in the following areas: fascial closures (including abdominal wall, thoracic and extremity fascial closures), gastrointestinal anastomoses, ophthalmic or microsurgery, cardiovascular, neural, osseus or tendinous tissue.

Healthcare professionals should be familiar with surgical procedures and techniques involving absorbable barbed device before employing STRATAFIX™ Spiral MONOCRYL™ Plus Device for wound closure, as risk of tissue separation / wound dehiscence leading to impaired healing may vary with the site of application and device material used. Healthcare professionals should consider the *in vivo* performance (under the Performance / Actions section) when selecting a device for use in patients.

Acceptable surgical practice should be followed for the management of contaminated or infected wounds.

As this is an absorbable device, the use of supplemental non-absorbable sutures should be considered by the healthcare professional in the closure of the sites that may undergo expansion, stretching or distension, or that may require additional support.

The use of STRATAFIX™ Spiral MONOCRYL™ Plus Device should not substitute normal observance of hygiene and/or otherwise needed antibiotic treatment. Healthcare professionals should consider the quantity and quality of tissue in which the loop will be anchored. The use of this device may be inappropriate in highly vascularized tissue or fragile tissue that cannot withstand potential further tightening/constriction within the loop.

Small bowel obstruction (SBO); including volvulus, bowel infarction, and significant morbidity, have been reported due to barbed suture ends hooking onto adjacent small bowel or omentery, such as in peritoneal closure.

Care should be undertaken to avoid leaving barbed suture ends adjacent to the peritoneum in extra-peritoneal tissue closure.

Avoid contacting STRATAFIX™ Spiral MONOCRYL™ Plus Device with other materials (e.g., surgical gauge, drapes, etc.) in the surgical field to prevent engraining on the bars. If the bars catch, carefully pull the material in the opposite direction of the needle to disengage it from the bars.

Do not resterilize/reuse. Reuse of this device (or portions of this device) may create a risk of product degradation that may result in device failure and/or cross-contamination that may lead to infection or transmission of blood-borne pathogens to patients and healthcare professionals.

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As an absorbable device, this device may act transiently as a foreign body. As may any foreign body, prolonged contact of any device with salt solutions, such as those found in the urinary or biliary tracts, may result in calculus formation. Like all foreign bodies, STRATAFIX™ Spiral MONOCRYL™ Plus Device may potentiate infection.

Discard unintentionally opened / partially used / used devices and packages.

Precautions

STRATAFIX™ Spiral MONOCRYL™ Plus Device contains a loop end and bars to anchor tissues and does not require knots to secure the device. Tying of knots on the barbed section of the device will damage the bars and potentially reduce their effectiveness and may result in extended surgery time or treatment failure and additional surgery. For the device to function properly, STRATAFIX™ Spiral MONOCRYL™ Plus Device must first be anchored by using the fixation loop in robust tissue. Then subsequently engage the bars into the tissue in a standard closure pattern. When completing placement of the barbed segment in subcuticular tissue, take at least one additional pass in the reverse direction to lock STRATAFIX™ Spiral MONOCRYL™ Plus Device in place.

Subcuticular devices should be placed as deeply as possible to minimize the erythema and induration normally associated with the absorption process.

Under some circumstances, notably orthopaedic procedures, immobilization of joints by external support may be employed at the discretion of the healthcare professional. Consideration should be taken in the use of absorbable devices in tissues with poor blood supply as device extrusion and delayed absorption may result.

This device may be inappropriate in elderly, malnourished or debilitated patients, or in patients suffering from conditions that may lead to wound healing.

In handling this or any other device material, care should be taken to avoid damage. Avoid crushing or crimping damage due to application of surgical instruments such as a forceps or needle holders. Do not attempt to remove memory in the polymer by running fingers down the device material as this can damage the bars.

Care should be taken to avoid tearing when handling surgical needles. Grasp the needle in an area one-third (1/3) to one-half (1/2) of the distance from the attachment end to the point. Grasping in the point area could impair the penetration performance and cause fracture of the needle. Grasping at the attachment end could cause bending or breakage. Reshaping needles may cause them to lose strength and be less resistant to bending and breaking.

Healthcare professionals should exercise caution when handling surgical needles to avoid inadvertent needle-stick injury that may result in transmission

of blood-borne pathogens from contaminated needles. Broken needles may result in extended or additional surgeries or residual foreign bodies. In the event of a product malfunction before use, such as a bent, broken or detached needle, or device damage including barbs, the product should be discarded and a new one obtained to begin the procedure. In the event of product malfunction during use, it is up to the discretion of the healthcare professional whether to continue or discontinue usage of the product and how to complete the procedure. Refer to Application / Instructions for Use section for instructions on removal of the device.

Adverse Reactions / Undesirable Side Effects

Adverse reactions associated with the use of this device include tissue separation / wound dehiscence leading to impaired healing or failure to provide adequate wound support in the closure of the sites where expansion, stretching, or distension occur, and in elderly, malnourished or debilitated patients or in patients suffering from other conditions that may delay wound healing. Other associated adverse events include calculus formation in urinary or biliary tracts when prolonged contact with salt solutions such as urine or bile occurs, transient inflammatory tissue reaction, allergic reaction (local hypersensitivity) to IRGACARE™ MP (triclosan), and delayed absorption of the wound segment near the needle, as the bars are required for successful wound approximation with STRATAFIX™ Spiral MONOCRYL™ Plus Device.

Device breakage may result in hemorrhage. Improper use of the device may cause soft tissue injury. Damage to device during use can cause extended surgery time or treatment failure and additional surgery.

Broken needles may result in extended or additional surgeries or residual foreign bodies. Inadvertent needle sticks with contaminated surgical needles may result in the transmission of blood-borne pathogens.

Healthcare professionals should convey adverse reactions, undesirable side effects and risks associated with the product and the procedure to the patient and address the risks and the procedure with the healthcare professional in case of any deviation from the normal postoperative course.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the country competent authority.

Magnetic Resonance Imaging (MRI) / Carcinogenic, Mutagenic, Toxic to Reproduction (CMR) / Endocrine Disrupting (ED) Safety Information STRATAFIX™ Spiral MONOCRYL™ Plus Device is MR safe. No known CMR Category 1a/1b and ED substances are present at >0.1%. Category 1a/1b are defined as known or presumed human carcinogens (H340), mutagen (H350) or reproductive toxicant (H360), based on human evidence and animal studies.

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